


# TIRED OF INJECTIONS?

Insulin Delivery For Kids Has  
Made A Huge DASH Forward



omnipod<sup>®</sup>  
DASH™

INSULIN MANAGEMENT SYSTEM



*Lexi B.*

PODDER™ SINCE 2013

**KID-FRIENDLY  
INSULIN  
THERAPY  
THAT'S  
SIMPLE,  
SMART,  
DISCREET**



*Trent W.*

PODDER™ SINCE 2013

Questions? Call us anytime  
at 800-591-3455 or visit  
[www.myomnipod.com](http://www.myomnipod.com)

**1 A Tubeless, Waterproof\* Pod**

Free your child from daily injections, tubing hassles, and limited clothing choices

**2 A Touch-Screen Personal Diabetes Manager (PDM)**

Discreet insulin delivery with a few finger taps

**3 iOS Mobile apps**

**Omnipod VIEW™** Allows caregivers to remotely view PDM data

**Omnipod DISPLAY™** Convenient access to Omnipod DASH™ System data from your child's smartphone


**Omnipod® iPhone widgets** A quick glance of PDM and CGM data together on one screen



Download on the App Store



Download on the App Store



Optimized For Use With  
**Contour  
next ONE**  
Blood Glucose  
Monitoring System

\* The Pod has an IP28 rating for up to 25 feet for 60 minutes. The PDM is not waterproof.

**SO SIMPLE TO  
GET STARTED**

*Taylor F.*

PODDER™ SINCE 2017

## **1** **FILL THE POD**

Fill Pod with up to 200 units of insulin

## **2** **APPLY THE POD**

The tubeless Pod can be worn almost anywhere an injection would be administered

## **3** **TAP “START” ON THE PDM**

The small, flexible cannula inserts automatically; they'll never see it and barely feel it



Questions? Call us anytime  
at 800-591-3455 or visit  
[www.myomnipod.com](http://www.myomnipod.com)

**DOSE  
EASILY AND  
DISCREETLY**

Questions? Call us anytime  
at 800-591-3455 or visit  
[www.myomnipod.com](http://www.myomnipod.com)



**Spend Less Time On Math,  
More Time Enjoying Meals**

Built-in bolus calculator

**Simplify Carb Counting**

Integrated CalorieKing® food library\* with  
nutrition facts for over 80,000 foods

**Adjust Basal Rates and Deliver  
Bolus Doses Quickly and Easily**

Familiar smartphone-like technology with  
simple-to-use display

**Multiple Languages**

Available in English and Spanish

\* CalorieKing® food library available in the English language only

# SMART MOBILE APPS MAKE IT EASY TO SEE YOUR CHILD'S DATA




## Omnipod VIEW™

Lets you and up to 12 loved ones view their insulin delivery information remotely



## Omnipod DISPLAY™

Your child can access their insulin delivery information on their smartphone and quickly find their PDM with the "Find My PDM" feature



## Omnipod® iPhone widgets\*

A combined view of your child's insulin delivery and CGM data on one screen



**Insulet Provided Glooko**

Gives healthcare providers better insight into your child's insulin management

- Eliminates data upload time during office visits
- Enables you and your healthcare providers to collaborate on treatment and make better diabetes management decisions

Questions? Call us anytime at 800-591-3455 or visit [www.myomnipod.com](http://www.myomnipod.com)


\* The Dexcom System does not have integrated functionality with the Omnipod DASH™ System. iPhone widgets are mobile application shortcuts that are visible after swiping right from the home screen to the iPhone's Today View

Omnipod DISPLAY™ and Omnipod VIEW™ apps are available with iOS 11 compatible smartphone devices iPhone 5S, 6, 6 Plus, 6S, 6S Plus, SE, 7, 7 Plus, 8, 8 Plus, and iPhone X.

Glooko is a trademark of Glooko, Inc. and used with permission.

# NO NEED FOR MANUAL BLOOD GLUCOSE ENTRIES

The Only Meter to Wirelessly Sync Readings to The Omnipod DASH™ Personal Diabetes Manager



*Aiden J.*

PODDER™ SINCE 2014

## OPTIMIZED FOR USE WITH THE OMNIPOD DASH™ SYSTEM, AND YOUR LIFE.

- The highly accurate<sup>1</sup> Bluetooth®-enabled CONTOUR®NEXT ONE Meter wirelessly transmits blood glucose readings to the Personal Diabetes Manager
- Eliminates the need for manual entry of blood glucose readings, and the potential for user entry errors when calculating insulin dosing
- The smartLIGHT target range indicator aligns with the PDM and immediately identifies if blood glucose readings are above, below or within the patient's target range<sup>2</sup>
- Second-Chance® sampling allows you 60 seconds to reapply more blood which may help prevent wasted test strips<sup>2</sup>
- Stores up to 800 blood glucose readings


Contour  
next ONE  
Blood Glucose  
Monitoring System



<sup>1</sup> Christiansen, M. et al. A New, Wireless-enabled Blood Glucose Monitoring System That Links to a Smart Mobile Device: Accuracy and User Performance Evaluation. Journal of Diabetes Science and Technology, 2017, Vol. 11(3) 567- 573.

<sup>2</sup> CONTOUR®NEXT ONE User Guide

# APPLY A NEW POD EVERY 3 DAYS AND AVOID 14 INJECTIONS



Average number of injections required every 3 days  
(MDI: ~4.5 per day)

## Consistent, Hands-Free Insertion

- No need for your child to see or touch an insertion needle
- They can wear the Pod almost anywhere their insulin would be injected manually

## Discreet Insulin Delivery

- Quick and easy bolus dosing delivered with a few finger taps on the Personal Diabetes Manager
- No need for your child to interrupt social activities for manual injections

## Precise

- Remove the guesswork and variability of injections
- Eliminate the need for blood glucose logs
- Deliver insulin in increments of 0.05 units



# FEEL WHAT IT'S LIKE TO WEAR A POD

Get a Podder™  
Experience Kit at  
No Cost – Just Call  
800-591-3455

# WHAT MAKES PODDERS™ HAPPIER THAN TRADITIONAL PUMPERS?

## Waterproof\*

- No need to disconnect to shower or swim
- Kids can stay connected to their insulin delivery and avoid disruption of their favorite activities

## Comfortable

- No need to make clothing choices to accommodate bulky pumps
- The Pod's auto-insertion feature allows for more convenient access to a variety of placement sites

## Tubeless

- No range of motion limits, ideal for active kids
- No concerns that tubes will pull out
- No struggles with infusion sets

## Zero Basal

- Children with low insulin requirements can set their basal segments as low as 0.00 units/hour

Questions? Call us anytime  
at 800-591-3455 or visit  
[www.myomnipod.com](http://www.myomnipod.com)



## ACTUAL SIZE



*Giuliana P.*

PODDER™ SINCE 2015

\* The Pod has an IP28 rating for up to 25 feet for 60 minutes. The Personal Diabetes Manager is not waterproof.

**FOR KATIE,  
LIFE AS A  
PODDER™  
IS SO  
MUCH MORE  
ENJOYABLE**

**SCHOOL,  
SPORTS,  
AND A LOT OF  
SILLINESS**

After six months of multiple daily injections, Katie loves how the Pod enables her to be more spontaneous. “We go to the lake in the summer so at any moment my friends and I could be swimming, tubing, kayaking. It’s just so much easier now.”

*“I don’t have time  
to slow down.  
With the Pod, I  
don’t have to.”*

*Katie R.*

PODDER™ SINCE 2015

Questions? Call us anytime  
at 800-591-3455 or visit  
[www.myomnipod.com](http://www.myomnipod.com)

# EXCEPTIONAL CUSTOMER SERVICE AND SUPPORT

## Product Support

Available to you 24 hours a day, 7 days a week, even if you are traveling abroad

## Multiple Languages

Translation services available for over 150 languages

## Patient Training

Access to Joslin ADEPT™ certified clinicians to support your educational needs\*

Find Out If The  
Omnipod DASH™ System  
is Covered By Your  
Insurance

Call **877-467-8538** or Visit  
**[omnipodinsurance.com](http://omnipodinsurance.com)**

Questions? Call us anytime  
at 800-591-3455 or visit  
**[www.myomnipod.com](http://www.myomnipod.com)**

\* Advanced Diabetes Education in Pump Therapy. ADEPT is a trademark of Joslin Diabetes Center and used with permission.



*Nathaniel P.*

PODDER™ SINCE 2014

# TECHNICAL SPECIFICATIONS

## Pod Design

### Size

1.53" wide x 2.05" long x 0.57" high

### Insulin

Holds up to 200 units of U-100 rapid acting insulin

### Weight

0.92 oz without insulin

## Personal Diabetes Manager Features

### Alarms

Adjustable volume and vibration alarms

### Alerts/Reminders

BG Reminder  
Pod Expiration  
Low Reservoir  
Bolus Reminders  
Program Reminders

## Other Features

**4.0" LCD Touch Screen**

**Bolus Calculator**

**CONTOUR®NEXT ONE  
Bluetooth® Wireless Technology**

**Holds Up to 90 Days of Data**

## Bolus Insulin Delivery

### Bolus Increments

0.05 units

### Correction Factor

1:1 to 1:400 mg/dL  
1/mg/dL increments

### Insulin to Carb Ratio

1:1 to 1:150 g with 0.1 g increments

### Maximum Bolus Size

30 units

## Basal Insulin Delivery

### Minimum Basal

0.00 units/hour

### Maximum Basal

30 units/hour

### Basal Increments

0.05 units/hour

### Temporary Basal Duration

30 minutes to 12 hours  
30 minute increments

## CONTOUR®NEXT ONE Blood Glucose Meter\*

### Size

97mm L x 28mm W x 14.9mm H

### Weight

36g

### Battery

Two 3v CR2032 or DL 2032 coin cell

### Battery Life

1000 tests (1 year average use;  
3 tests per day)

### Memory

Stores 800 most recent readings

### Test Strips

CONTOUR®NEXT Test Strips

### Test Sample

Capillary whole blood only

### Sample Volume


0.6µL

\* Included in the Omnipod DASH™ System Welcome Kit



*Frankie C.*

PODDER™ SINCE BEFORE  
HE COULD EVEN WALK



omnipod®  
DASH™

INSULIN MANAGEMENT SYSTEM

Customer Care: 800-591-3455

[myomnipod.com](https://myomnipod.com)

Insulet Corporation  
600 Technology Drive  
Suite 200  
Billerica, MA 01821

© 2018 Insulet Corporation. Omnipod, the Omnipod logo, DASH, the DASH logo, Omnipod DISPLAY, Omnipod VIEW and Podder are trademarks or registered trademarks of Insulet Corporation. All rights reserved. Ascensia, the Ascensia Diabetes Care logo, and Contour are trademarks and/or registered trademarks of Ascensia Diabetes Care. Glooko is a trademark of Glooko, Inc. and used with permission. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Insulet Corporation is under license. Dexcom and Dexcom G6 are registered trademarks of Dexcom, Inc. and used with permission. All other trademarks are the property of their respective owners. The use of third party trademarks does not constitute an endorsement or imply a relationship or other affiliation.

40214-ENG AW R4 10/18